

FGF2 alternative peptide (FGFR1c agonist)

Code	PG-011
MW (Acetate):	5127.81
Purity:	≥95% (HPLC)
Amount:	10 µg

Formulation: This vial contains 10 µg of the titled compound, which has been lyophilized from 0.2 µm filtered DMSO solution.


Reconstitution: Reconstitute at 100~200 µg/mL in **DMSO**. Centrifuge the vial before opening.
Rinse the inside of the vial thoroughly with DMSO.
 When further diluting the DMSO solution, we recommend using a buffered solution (e.g., PBS) with a carrier component to prevent peptide adhesion.
 <Example>
 A 100 µg/mL solution can be prepared by reconstituting 10 µg of peptide in 100 µL DMSO.

Storage: Store at -20°C or below for the container vial unopened. Recommend to use the peptide as quickly as possible after opening and reconstitution.


Activity: Measured by its ability to induce phosphorylation of human FGFR1c in BaF3 cells overexpressing hFGFR1c in the presence of 1 U/mL of heparin, as well as through a proliferation assay using human mesenchymal stem cells (MSCs) derived from bone marrow, conducted without heparin.
The expected ED50 for these effects is 0.1–1.5 ng/mL.

Evaluation Guideline: PG-011 is expected to **exhibit comparable activity at 1/3(ng/mL basis) of rFGF2 concentration.** For MSC or fibroblast proliferation, **test at several points around 1/3 of rFGF2.** For iPSC maintenance, **start at 1/3 and test around higher concentrations.** For differentiation induction, **start at 1/3 and test around lower concentrations.** Optimize by testing multiple concentrations for each system.

Phosphorylation of FGFR1c



Proliferation of human MSCs



This product is "For Research Use Only"